# **FX-i/FZ-i**

**Quick Start Guide – FX-i/FZ-i Series  
Compact Precision Balances**

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

 This is a hazard alert mark

## UNPACKING

- The balance is a precision instrument. Unpack the balance carefully. Keep the packing material to be used for transporting the balance in the future.
- The packing contents depend on the balance model. See the illustrations to confirm that everything is contained.

 Please confirm that the AC adapter type is correct for your local voltage and receptacle type.

Quick Start Guide

# SETTING UP

### Setting up your balance


1. Place the balance on a solid and level weighing table.
2. The FZ-i series (all) and FX-120*i* /200*i* /300*i* /500*i*:  
Assemble the pan support, weighing pan and breeze break, on the balance as shown in the illustration.  
  
FX-1200*i* /2000*i* /3000*i* /5000*i*:  
Assemble the pan support and weighing pan on the balance as shown.
3. Adjust the leveling feet to level the balance. Confirm it using the bubble spirit level.
4. Confirm that the AC adapter type is correct for the local voltage and power receptacle type.
5. Connect the AC adapter to the balance. Warm up the balance for at least 30 minutes with nothing on the weighing pan.

**Note:**

*The AC adapter plug is protected against dust and may be difficult to insert. When inserting the plug, turn the plug while pushing on it.*

### Balance location

- Install the balance in an environment where the temperature and humidity are not excessive. The best operating temperature is about 20°C / 68°F at about 50% relative humidity.
- Install the balance where it is not exposed to direct sunlight and it is not affected by heaters or air conditioners.
- Install the balance where it is free of dust.
- Install the balance away from equipment which produces magnetic fields.
- Install the balance in a stable place avoiding vibration and shock. Corners of rooms on the first floor are best, as they are less prone to vibration.
- The weighing table should be solid and free from vibration and drafts, and as level as possible.
- Level the balance by adjusting the leveling feet and confirm it using the bubble spirit level.
- Ensure a stable power source when using the AC adapter.
- Plug in the AC adapter and warm up the balance for at least 30 minutes.
- Calibrate the balance before use or after having moved it to another location.

### Caution

**Do not install the balance where flammable or corrosive gas is present.**

---


## PART NAMES AND FUNCTIONS

---


### Key operation

Key operation affects how the balance functions. The basic key operations are:

- “Press and release the key immediately” or “Press the key”  
= normal key operation during measurement
- “Press and hold the key”

### Display symbols

Each key, when pressed or when pressed and held, functions as follows:

Key	When pressed	When pressed and held
	Turns the display on or off. The standby indicator is displayed when the display is turned off. The weighing mode is enabled when the display is turned on. This key is available anytime. Pressing the key during operation will interrupt the operation and turn the display off.	
	In the weighing mode, turns the minimum weighing value on or off. In the counting or percent mode, enters the sample storing mode.	Enters the function table mode. Refer to “10. FUNCTION TABLE” in the full instruction manual.
	Switches the weighing units stored in the function table. Refer to “5. WEIGHING UNITS” in the full instruction manual.	Enters the response adjustment mode.
	Cancels the operation when performing settings. With the FZ-i series, enters the calibration mode using the internal mass. (One-touch calibration.)	Enters the calibration mode.
	Outputs the weight data to a printer or personal computer using the RS-232C serial interface, depending on the function table settings. Confirms the operation when performing settings.	No function at the factory setting By changing the function table: Outputs “Title block” and “End block” for GLP report.
	Sets the display to zero.	

---

# CALIBRATION

---

### Calibration Using the Internal mass (One-Touch Calibration FZ-i)

This function calibrates the balance using the internal mass. The only operation required is to press the **CAL** key.

#### Operation

1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the weighing pan.
2. Place the accessory breeze break on the balance and press the **CAL** key.
3. The balance displays **Calin** and performs calibration using the internal mass. Do not allow vibration or drafts to affect the balance.
4. The balance displays **end** after calibration. If the "GLP output (info)" parameter of the function table is set to "1" or "2", the balance displays **gfp** and outputs "Calibration Report" using the RS-232C interface or stores the data in memory. For details on the calibration report format, refer to "11-2 GLP Report".
5. The balance will automatically return to the weighing mode after calibration.

#### About the Internal Mass

The mass of the internal mass may change due to corrosion or other damage caused by the operating environment, or due to aging. Check the internal mass periodically. Correct the internal mass value as necessary. For details, refer to the user manual.

To maintain the weighing accuracy, perform the internal mass value correction periodically.

## Quick Start Guide

### Calibration Using an External Weight

This function calibrates the balance using an external weight.

#### Operation

1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the pan.
2. Press and hold the **CAL** key until **Calout** is displayed, and then release the key.
3. The balance displays **Cal 0**.

- If you want to change the calibration weight, press the **SAMPLE** key and go to step 4.  
➤ If you use the calibration weight value stored in the balance, go to step 5.

4. Specify the calibration weight value as follows:

<b>SAMPLE</b> key	To toggle the display between calibration weight selection mode or calibration weight value adjustment mode.
<b>RE-ZERO</b> key	(To increase the value.)
<b>MODE</b> key	(To decrease the value.)
	To select the calibration weight or adjust the value.
<b>PRINT</b> key	To store the new weight value. Even if the AC adapter is removed, the data is maintained in non-volatile memory.
<b>CAL</b> key	To cancel the operation and return to <b>Cal 0</b> .

5. Confirm that there is nothing on the pan and press the **PRINT** key. The balance measures the zero point. Do not allow vibration or drafts to affect the balance.

The balance displays the calibration weight value.

6. Place the displayed calibration weight on the pan and press the **PRINT** key. The balance measures the calibration weight. Do not allow vibration or drafts to affect the balance.

7. The balance displays **end**. Remove the weight from the pan.

8. When the GLP output is set, the balance displays **glp** and outputs "Calibration Report".

9. The balance will automatically return to the weighing mode.

10. Place the calibration weight on the pan and confirm that the value displayed is within  $\pm 2$  digits of the specified value. If it is not within the range, check the ambient conditions such as drafts and vibration, also check the weighing pan. Then, repeat steps 1 to 9.

---

## ERROR CODES

---

### Error Codes

Display	Error code	Description
	EC, E11	<p><b>Stability error</b>            The balance cannot stabilize due to an environmental problem. Prevent vibration, drafts, temperature changes, static electricity and magnetic fields.            Refer to "3. PRECAUTIONS" in the full instruction manual for details on the operating environment and "7. RESPONSE ADJUSTMENT" in the full instruction manual about adapting the balance to the environment.            To return to the weighing mode, press the <b>CAL</b> key.</p>
		<p><b>Out of range error</b>            The value entered is beyond the settable range.            Re-enter the value.</p>
	EC, E17	<p><b>Internal mass error (Only for the FZ-i series)</b>            The internal mass application mechanism does not function properly. Perform the weighing operation from the beginning again.</p>
	EC, E20	<p><b>Calibration weight error</b>            The calibration weight is too heavy.            Confirm that the weighing pan is properly installed.            Confirm the calibration weight value.            Press the <b>CAL</b> key to return to the weighing mode.</p>
	EC, E21	<p><b>Calibration weight error</b>            The calibration weight is too light.            Confirm that the weighing pan is properly installed.            Confirm the calibration weight value.            Press the <b>CAL</b> key to return to the weighing mode.</p>
		<p><b>Overload error</b>            A sample beyond the balance weighing capacity has been placed on the pan.            Remove the sample from the pan.</p>
		<p><b>Weighing pan error</b>            The weight value is too light.            Confirm that the weighing pan is properly installed and calibrate the balance.</p>
		<p><b>Sample mass error</b>            The balance cannot store the sample for the counting mode or for the percent mode because it is too light.            Use a sample that is heavier.</p>

# Quick Start Guide

# FX-i/FZ-i Series Compact Precision Balances

Display	Error code	Description
  		<p><b>Unit mass error</b>            The sample unit mass for the counting mode is too light. Storing and using it for counting will cause a counting error.            Add samples to reach the specified number and press the <b>PRINT</b> key.</p> <p>Pressing the <b>PRINT</b> key without adding samples will shift the balance to the counting mode. But, to acquire accurate weighing, be sure to add samples.</p>
		<p><b>Balance internal error</b>            If this error appears persistently, contact the local A&amp;D dealer.</p>
	EC, E00	<p><b>Communications error</b>            A protocol error occurred in communications.            Confirm the format, baud rate and parity.</p>
	EC, E01	<p><b>Undefined command error</b>            An undefined command was received.            Confirm the command.</p>
	EC, E02	<p><b>Not ready</b>            A received command cannot be processed.            e.g. The balance received a Q command, but not in the weighing mode.            e.g. The balance received a Q command while processing a RE-ZERO command.            Adjust the delay time to transmit a command.</p>
	EC, E03	<p><b>Timeout error</b>            If the timeout parameter is set to "t-Up1", the balance did not receive the next character of a command within the time limit of one second.            Confirm the communication.</p>
	EC, E04	<p><b>Excess characters error</b>            The balance received excessive characters in a command.            Confirm the command.</p>
	EC, E06	<p><b>Format error</b>            A command includes incorrect data.            e.g. The data is numerically incorrect.            Confirm the command.</p>
	EC, E07	<p><b>Parameter setting error</b>            The received data exceeds the range that the balance can accept.            Confirm the parameter range of the command.</p>
<b>Other errors</b>		If the errors described above cannot be released or other errors are displayed, contact the local A&D dealer.

### Legal-for-Trade Model Information

#### NOTE FOR EC APPROVAL MODELS

FZ-120i-EC, FZ-200i-EC, FZ-300i-EC, FZ-1200i-EC, FZ-2000i-EC, FZ-3000i-EC

Model	FZ-120i-EC	FZ-200i-EC	FZ-300i-EC	FZ-1200i-EC	FZ-2000i-EC	FZ-3000i-EC
Class	II	II	II	II	II	II
Max	122g	220g	320g	1220g	2200g	3200g
Min	0.02g	0.02g	0.02g	0.5g	0.5g	0.5g
Min Display d	0.001g	0.001g	0.001g	0.01g	0.01g	0.01g
Min Display e	0.01g	0.01g	0.01g	0.1g	0.1g	0.1g

Class	Item and Parameter			Description
Func	trc Zero tracking	0	OFF	Keeps zero display by tracking zero drift.
		1	Normal	
		2	Not Available	
		3	Not Available	
dout Data output	ar-d Zero after output		Not Available	Adjusts zero automatically after data is output
C5 in 1 *	Internal mass value correction: Method 1			Displayed only when the function switch (internal mass value correction switch) is set to 1.
C5 in 2 *				
Internal mass value correction: Method 2				

- Calibration using an external weight (CalOut) is not available.
- Automatic span adjustment will be enabled (the balance will automatically calibrate using the internal weight due to environmental temperature changes).

**Note about Units:** When EU 2014/31/EU verification is performed only the following units will be available.

Unit Name	Unit Symbol	Unit Name	Unit Symbol
Grams / Milligram	g / mg	Pieces	pcs
Programmable Unit	MLT	Percent	%

## Notes about Verification

### Metrological note for EU countries

The identification plate of the scale has a CE mark followed by the metrology mark (M + year two-digit; framed) and the number of the notified body. Such a scale can be started up and operated on-site provided that it is not connected to an auxiliary device which is subject to metrological approval. Scales which are connected on site to an auxiliary device which is subject to metrological approval must undergo a conformity assessment or metrological approval. Only then, legal-for-trade use is permissible.

Scales without metrology mark or control marks may not be used in legal-for-trade applications.

Scales with an accuracy – Class II, with 'n' greater than 3300 have been adjusted to the place of use and may not be used elsewhere.

As per legal requirements, the scale user must ensure that the scale is used for its intended purpose. This includes the observance of legal requirements, especially in case of changes, amendments and extensions.

The verification of a scale is executed according to the respective statutory country regulations. The verification validity starts as soon as the device is put on the market.

---

## DISPOSAL OF THE DEVICE

---

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. To ensure that the device is disposed of correctly please recycle in accordance with local regulations. For further information contact [REDACTED] the supplier from whom you purchased the product.

---

## **EU Declaration of Conformity**

---

This is a Declaration of Conformity for the instrument below.

INSTRUMENT MODEL	FZ-i-EC Series
NAME AND ADDRESS OF THE MANUFACTURER	<b>A&amp;D Instruments Ltd</b> 24-26 Blacklands Way Abingdon Business Park Abingdon Oxfordshire OX14 1DY
TYPE APPROVAL CERTIFICATE	T7678
NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE B: EU TYPE EXAMINATION (Annex II .1 of 2014/31)	NMI 0122
NOTIFIED BODY ISSUING THE APPROVAL FOR MODULE D: CONFORMITY TO TYPE BASED ON QUALITY ASSURANCE OF THE PRODUCTION PROCESS (Annex II.2 of 2014/31)	SGS UK 0120

Conformity to the following Directives is demonstrated by reference to the following harmonized standards or technical documentation:

APPLICABLE DIRECTIVES	HARMONIZED STANDARDS OR OTHER TECHNICAL DOCUMENTATION
Directive on non-automatic weighing instrument (2014/31/EU)	EN 45501:2015
Directive on Electromagnetic Compatibility (2014/30/EU)	EN 61326-1:2013 Electrical equipment for measurement control and laboratory use – EMC requirements Part 1: General Requirements
Directive on Electrical Equipment designed for use within certain voltage limits (2014/35/EU)	EN 60950-1:2006 +A11:2009+A1:2010+A12:2011+A2:2013 Safety of Information Technology Equipment
Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II) (2011/65/EU)	EN 50581:2012

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

Signed for and on behalf of:	<b>A&amp;D Instruments Ltd</b>
Place and date of issue:	Abingdon, UK, 1 <sup>st</sup> October 2018
Name, function:	J. Ghuman, Managing Director
Signature:	

---

## **EU Declaration of Conformity**

---

**A&D Instruments Ltd. hereby declare that the following Weighing product conforms to the requirements of the council directives on ...**

Electromagnetic Compatibility (EMC) 2014/30/EU,

Low Voltage Equipment (LVD) 2014/35/EU and

Restriction of the use of certain Hazardous Substances (RoHS) 2011/65/EU

**provided that they bear the CE mark of conformity.**

**Model/Series....FX-i / FZ-i Series**

Standards applicable:

**EN 61326-1:2013**

Electrical equipment for measurement, control and laboratory use -EMC requirements Part 1: General requirements

**EN 60950-1:2006/A11:2009/A1:2010/A12:2011/A2:2013**

Safety of Information Technology Equipment

**EN 50581:2012**

Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

This Declaration of Conformity is issued under the sole responsibility of the manufacturer:

A&D Instruments, 24-26 Blacklands Way, Abingdon, Oxfordshire, OX14 1DY, England

**CE Mark first applied 08 July 2002**

Signed for A&D Instruments on 1<sup>st</sup> October 2018

A handwritten signature in black ink, appearing to read 'J. Ghuman'.

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