

**Supported OS Versions for Windows Communication Tools Software** (as of February 05, 2018)

Free Software						
		Balance/Scale		Multiple Balances and/or Scales	AD-4212C Series	
		WinCT Ver.2.01 (CD-ROM)	WinCT Ver.3.00	WinCT Ver.5.40	WinCT-Plus Ver.3.00	WinCT-AD4212C Ver.1.02
		RsCom/Key	RsCom/Key/Weight	RsCom/Key/Weight	RsMulti	
Win 98 to ME		Yes	Yes	No	Yes	Yes
Win XP		Yes	Yes	No	Yes	Yes
Win Vista		No	No	Yes	Yes	Yes
Win 7	32bit	No	No	Yes	Yes	Yes
	64bit	No	No	Yes	Yes	Yes
Win 8.1	32bit	No	No			
	64bit	No	No	Yes	Yes	Yes
Win 10	32bit	No	No			
	64bit	No	No	Yes	Yes	Yes
Remarks	CD-ROM included with the product as standard (discontinued in Apr 2014)		Available for download from the support page on the A&D website		Comes with the AD-8526 / Download from the A&D website * Microsoft Excel cannot be started from RsMulti.	

Paid Software				Serial/USB Converters	
		Pipette Accuracy Tester	Moisture Analyzer	Viscometer/Rheometer	AX-USB-9P/AX-USB-25P/AX-USB-DIN
		WinCT-Pipette Ver.1.26	WinCT-Moisture Ver.2.41	WinCT-Viscosity Ver.1.13	Device Driver *1
			RsFig/Temp	RsVisco	
Win 98 to ME		Yes	Yes	Yes	Yes *2
Win XP		Yes	Yes	Yes	Yes
Win Vista		Yes	Yes	Yes	Yes
Win 7	32bit	Yes	Yes	Yes	Yes
	64bit	Yes	Yes	Yes	Yes
Win 8.1	32bit				Yes
	64bit	Yes	Yes	Yes *	Yes
Win 10	32bit				Yes
	64bit	Yes	Yes	Yes	Yes
Remarks			* Some characters are not displayed.		*1 For Win 7/8.1/10, download from the ATEN website at <a href="http://www.aten.com/download/download_dv.php">http://www.aten.com/download/download_dv.php</a> (select the model No. UC232A) *2 Only supported with Win 98 SE.

Not verified yet

Note 1) WinCT-Moisture/Viscosity RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balance/Scale.

Note 2) WinCT-Pipette is included with Pipette Accuracy Testers.