

## Windows 10<sup>\*1</sup> Supported Versions of Communication Tools Software (as of February 05, 2018)

\*1 Verified with 64-bit versions only.

Free Software					
	Balance/Scale			Multiple Balances and/or	AD-4212C Series
	WinCT			WinCT-Plus	WinCT-AD4212C
	RsCom	RsKey	RsWeight	RsMulti	
Ver.	<b>5.40</b>	<b>5.40</b>	<b>5.40</b>	<b>3.00</b>	<b>1.02</b>
Win10 1511	Yes	Yes	Yes	Yes	Yes
Win10 1703	Yes	Yes	Yes	Yes	Yes
Win10 1709	Yes	Yes	Yes	Yes	Yes

Paid Software								
	Moisture Analyzer				Viscometer/Rheometer			Pipette Accuracy
	WinCT-Moisture <sup>*2</sup>				WinCT-Viscosity			WinCT-Pipette <sup>*4</sup>
	RsFig	RsTemp	RsCom	RsKey	RsVisco	RsCom	RsKey	
Ver.	<b>2.41</b>	<b>2.31</b>	<b>2.50</b>	<b>2.50</b>	<b>1.13</b>	<b>2.50</b>	<b>2.50</b>	<b>1.26</b>
Win10 1511	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1703	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1709	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes

\*2 Provided as standard for the MX-50 and MS-70

\*3 WinCT-Moisture/Viscosity RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balance/Scale.

\*4 Included with Pipette Accuracy Testers